

BEVERAGE CARE

SCIENTIFICALLY PROVEN Trusted By Beverage Professionals

MILKLINE POWDER

Performance Focus

Industry-leading research and development, combined with extensive testing ensures that our MILKLINE POWDER can be used across a whole spectrum of machines. This product is designed to be used manually or through a machine's automatic cleaning system to prevent component blockage by eliminating milk build up in pipe work and machine components to ensure a hygienically safe machine.

Our intensive focus on performance delivers a unique product providing the following significant benefits:

- SUPERIOR CLEANING ACTION Our MILKLINE POWDER contains high levels of active agents ensuring a superior cleaning action and a hygienically safe machine, while reducing cleaning times.
- 0 NON DAMAGING Product will not damage or degrade machine components unlike the more traditional hazardous products.
- ENVIRONMENTALLY FRIENDLY Proven by independent environmental audits to be non-toxic, biodegradable and non accumulative in the eco system.

Where to use!

- o BEAN TO CUP MACHINES
- o SERVICE ENGINEERS & WORKSHOP

When to use!

Use daily in accordance with machine manufacturers cleaning procedure.

How to use!

AUTOMATIC CLEANING

Pour 2-4 scoops into milk container and following machine manufacturer cleaning instructions.

MANUAL CLEANING

Pour 2-4 scoops into jug with 1 litre of hot water put milk pipe into jug and follow machine manufactures cleaning instructions.

Client Testimonial

SERVICE ENGINEER

'This product is super quick at dissolving and does not leave any undissolved powder when mixed with warm water. This means it is quick to use and I can be confident that no bits of powder remain undissolved leading to the risk of machine blockages when I leave site.'

Available Products

DISPENSING OPTIONPACK SIZECODE1 KILO TUB6 UNITS058354

Contact Bioguard today and experience the benefits we can bring to your operation

bioguardhygiene.com

tel +44 (0) 1604 592048 - email info@bioguardhygiene.com 64a St James Mill Rd, St James, Northampton, NN5 5JP, UK